



ALLOGRAFT SHIPPING & STORAGE SOLUTIONS




Suspension Frame & Insulated Shipper

Our suspension frame offers ultimate protection during the shipment and storage of allografts. Enveloped tightly in a durable film that suspends products within a frame, the suspension film absorbs shocks during impact, allowing allografts to move within a cushion of air. This prevents damage by eliminating contact with the insulated shipper.

- Reusable and single-use versions
- Recyclable
- Fully customizable






Vertical Insulated Shipper

Constructed with recyclable components and available in both reusable and single-use versions, our insulated shipper takes weight off of sensitive allograft implants by suspending them vertically inside of our customizable frame.

- Compatible with dry ice, wet ice and phase change materials
- Up to 120 hours of temperature control
- Reduces shock by 60%



LN₂ Freezer Rack Expansion Sleeves

Our customizable freezer rack expansion sleeves can double the capacity of your liquid nitrogen freezer.

- Labeled
- Recyclable
- Constructed with particulate free materials
- Fully customizable
- Remains flexible inside freezer



UFP
TECHNOLOGIES

100 Hale Street
Newburyport, MA 01950

800 372 3172
ufpt.com

Shaping Innovation

REV. 7.2017